

Home Phototherapy Best Practices

Wear UV Blocking Eye Protection.

Remove Skin Creams Prior To Treatment.


Cover or use SPF 30+ on unaffected body areas.

3 treatments per week is typical. Consult your physician.

Treatment Positioning Guide:

Consistent body positioning is key to a successful treatment. Use the distances outlined below:

200/400/600/600 LT/ 600 3D



200/400/HF/100/100 XL



Doses can be manually adjusted downward, but never below your starting dose. Selecting "Yes, 24-48 hours" will repeat the last dose with no increase.



To prevent overexposure, Guided mode has a lockout time of 16 hours before allowing the next treatment session.

Security Pin: Used to access the device and will always be 0007.

C-Code: A 4-digit code provided by your doctor, beginning with the letter C, used to refill exposures on devices that have exposure limiting software.


Getting Started With Guided Mode

STEP 1: Security PIN Entry

- Awaken the device by tapping the screen
- Tap the logo screen to continue
- Press "7" to display 0007, followed by the RIGHT ARROW 




SEE THIS SCREEN? CONTACT YOUR HEALTHCARE PROVIDER



Contact your doctor with the **C Code** shown on this screen. Your doctor will then provide you with a 4-digit exposure refill code. Enter the code followed by the RIGHT ARROW 

STEP 2: Confirm Prescription Details

- Verify the screen states the correct details
- Press the RIGHT ARROW  to proceed



This screen shows your prescribed starting dose. Press the RIGHT ARROW  to continue.

Future treatment doses will be adjusted based upon your feedback to these questions.



STEP 3: Prepare for Treatment

- Position your body – refer to the Helpful Tips section for correct positioning
- Put on your UV eye blocking protection
- Cover or use SPF 30+ on unaffected body areas.
- Press the RIGHT ARROW  to begin

