



STAND ASSEMBLY USER MANUAL

Tools Needed: Phillips-head screwdriver (“+” shaped end)

Phothera 200 Stand Accessory Kit Includes: 1 Cradle, 2 Spacer Disks, 2 Knobs, 4 Phillips Screws, and 2 Plastic Washers

1. Unplug the power cord from the device.
2. Find a suitable workspace and stand the device upright so that the lamps are facing toward you, and the timer is towards the top, shown in the illustration to the right.
3. Place one of the round spacer disks on the top of the device and align the two holes on the disk with the two holes in the device exterior.
4. Insert two screws into the holes and tighten with a Phillips screwdriver.
5. Place the other round spacer disk on the bottom side of the device and align with the holes there.
6. Insert two screws into the holes and tighten with a Phillips screwdriver.
7. Place the cradle next to the device and align the holes on the cradle with the holes on the spacer disks.
8. Insert the black knobs into the plastic washers and cradle holes. Screw in to tighten. Your device is now ready for use, as shown in the illustration on the right.

