

Phoothera

600 3D



Phoothera 600-3D

The Phoothera 600-3D is a full-body phototherapy device that encompasses the patient in 180 degrees of phototherapeutic light. The device offers a slim design and minimal footprint, variable lamp configurations, and can be equipped with either the DT or ClearLink Controller.

Phothera

600 3D



- 1 The Phothera 600-3D offers patients and prescribers choices when choosing lamp quantities and wavelengths:
 - 8 or 12 lamps
 - Narrowband UVB or UVA
- 2 The slim lines of the Phothera 600-3D make it ideal for a home setting or offices with minimal space; easily fitting within just 3 square feet of floor space.
- 3 To enhance efficiency and shorten treatment times, the doors on the 600-3D wrap the patient in 180 degrees of light. When closed, the doors discreetly cover the fluorescent lamps.
- 4 The Phothera 600-3D offers both simplicity and flexibility, with two customizable controller options designed to fit each patient's unique needs and goals.
 - ClearLink Controller - Safety, security, and patient monitoring, all in an intuitive, touch-screen controller.
 - DT Controller - Simplicity in the form of a time-only controller.

Treatment Area	33" x 70.25"	83.8 cm x 178.4 cm
Width (with doors open)	39.5"	100.3 cm
Height	74.5"	189.2 cm
Depth (including feet)	24"	61 cm
Weight	129 lbs	58.5 kg

For further information, please contact Phothera:

216.831.0600

phothera.com | info@phothera.com